

Veterinary Doppler Blood Pressure Monitor

Model: DopMeter 100

Features:

- 1) High quality manometer
- 1) Ultrasonic probe in 8.0MHz
- 2) Five pieces of oxford fabric cuffs
- 3) Ultrasonic host
- 4) Ultrasonic gel
- 5) Earphone
- 6) Power adapter
- 7) Packed in aluminum alloy box



Specifications:

Measuring range: 0mGs-199.9mGs
 Frequency of the probe: 8.0MHz
 Size of the blood pressure cuff: 4.2cm-7.1cm, 5.4cm-9.1cm, 6.9cm-11.7cm, 9cm-13cm, 12cm-16cm
 Rechargeable battery voltage: 7.4V, 1500mAh
 AC input: 100-240V, 50/60Hz
 DC output: DC8.4V, 1000mA
 Continuous working time: > 8 hours
 Speaker: 4Ω, 1W

Veterinary Doppler Blood Pressure Monitor

Model: DopMeter 200D/200C

Features:

- 1) Aneroid Sphygmomanometer, no pollution
- 2) Integrated sphygmomanometer, easy to operate
- 3) Specialized blood flow probe with high sensitivity
- 4) Abundant vascular doppler menu
- 5) Multiple size cuff, feasible to different animals
- 6) Display vascular sound volume and parameter
- 7) Turn off automatically when no using in 3 minutes
- 8) Both of battery and AC power supply

Applications:

Test the animal blood pressure in animal hospitals or at home.



Specifications:

Blood pressure range: 0-300mmHg
 Pressure precision: ±3.75mmHg
 Vascular probe frequency: 8Mhz
 Power supply: Built-in rechargeable battery, AC 220-110V, 50/60Hz

Standard:

Vascular doppler, 8MHz Vascular probe, Sphygmomanometer, Cuff, Adaptor, Carrying case.

Optional:

Welch Allyn sphygmomanometer.

Note:

DopMeter 200D: Black/white display screen.
 DopMeter 200C: Color touch screen.

