

Veterinary Micro Bone Saw

Model: DS-1411

Features:

- 1) Autoclavable to 135°C with 1.6mpa
- 2) Whole power is steady with low noise in operation
- 3) Handpiece can drive a variety quick connects
- 4) Attachments can be replaced easily
- 5) 360° Rotating
- 6) Main parts made of high-strength stainless steel

**Specifications:**

Hand piece: Working Voltage: 7.2V; Output Power: $\geq 50W$; Unit rise in temperature: $\leq 20^{\circ}C$;
Charger: Input Voltage: AC 110V-240V, 50Hz-60Hz; Output Power: 7.2V; Input Power: 30VA;
Battery: Voltage: 7.2V; Charging time: About 3 hours;
Saw: Frequency: 0-20000 os/m; Non-load noise: $\leq 45dB$; Output power: 50W; Sterilized: $\leq 135^{\circ}C$.

Standard Configuration:

1*Hand piece, 2*Saw blades, 2*Batteries, 1*Charger, 2*Sterilizing channel, 1*Aluminium case.