

Veterinary Otoscope

Model: ENT-EI

Application:

Otoscope is applied to ear diagnosis and general surgery.

Specifications:

Size of the aural: diameter 3mm, 4mm, 5mm, 6mm, 8mm
Magnifier: 3 times
Power Supply: 1.5Vx2 (C Type)
Illumination Bulb: 2.5V/0.3A
Battery handle: Plastic



Veterinary Direct Ophthalmoscope

Model: ENTA-I

Application:

Direct ophthalmoscope is applied for eye check.

Specifications:

Illumination form: big spot, small spot, red free
Diopter Compensation: 0 \pm 1 \pm 2 \pm 3 \pm 4 \pm 5 \pm 6 \pm 8 \pm 10 \pm 12 \pm 15 \pm 20 -25-30(D)
Power Supply: 1.5Vx2 (C Type)
Illumination Bulb: 2.5V/0.3A



Animal Anesthetic Laryngoscope

Model: AL-5

Specifications:

Length of lens:
No. 0: 75mm
No. 1: 100mm
No. 2: 155mm
No. 3: 195mm
No. 4: 250mm
Light source: LED/Halogen
Material: Stainless Steel
Power: two 1.5V batteries
Weight: 1.4kg

Optional:

AL-5A: LED lamp
AL-5B: Halogen lamp



Digital HD Veterinary Otoscope

Model: EndoBox-EarView

Features:

Multi-functional Integration (Imaging, Camera, Light Source, Record)
Universal compatible with optical lens
Small volume, Light weight, Portable
Support endoscopic mobile application
Various display, Various application

Specifications:

HD Monitor: 12 inches, 16:9
HD Camera: 1/2" Full HD CMOS Sensor
Video Output: HDMI, Synchronous
LED Light Source: Power 30W, Illuminance 700,000Lux
HD Record: Resolution 1920x1080, Format JPEG
Functions: Freeze, Replay, Snapshot, Storage
System Control: Button, Remote

Standard:

Digital HD Veterinary Endoscopy Camera System (12-inch Main unit EndoBox-12, Camera head, Optical adapter, Light cable, Remote-control, Power adapter, Power cable), Veterinary otoscope (C0010: 0°, Φ 2.7x107.5mm) and Veterinary Sheath (C1001).

Optional:

15 / 22 / 26 / 32-inch screen, Carrying case, Stand

Functions:

HD imaging, HD camera, Powerful LED light source
HD recording, Storage, Playback, Snap
Synchronous secondary output
Real-time network output, Rebroadcast
Network real-time browsing

Extensible Applications:

Laparoscope, Cystoscope, Arthroscope, Sinuscope, Laryngoscope, Thoracoscope, Flexible Fiberscopes

