

Animal Defibrillator

Model: DM8A



Feature:

- 1) Non-synchronizer: Outlife difibrillator
- 2) Sine wave: Monophasic technology
- 3) Defibrillator type: Manual, Asynchronized
- 4) Energy accuracy: <+1%
- 5) Energy selection: 0, 3, 5, 7, 10, 20, 30, 50, 100, 200, 300, 360 joules (nominal at a resistance of 50)
- 6) Charging time: Maximum 10 seconds to 360 Joules
- 7) Standard Adult/Pediatric Paddles Reusable external adult paddles (pediatric paddles integrated)
- 8) Power requirements:
 - a) AC Power supply 100V~240V,50Hz/60 Hz
 - b) Vehicle Voltage (Optional) DC 12V



Optional:

- Biphasic (DM8A-II)
- Li-ion battery
- Vehicle voltage

