

12.1" Animal Monitor

Model: Moni 6DV



Feature:

- 1) Colorful 12.1" TFT screen, backlight buttons available
- 2) Flexible, light, easy operating
- 3) Various language operating systems
- 4) Built-in detachable rechargeable lithium battery, 2.5 hours battery life
- 5) Supporting for wireless or wired networking

Alarming

- 3 level light/voice alarm, dual alarm light for physiological and technical alarming

Data Storage

- 1) 400 groups NIBP measurements
- 2) 72 hours trend data and trend graphs
- 3) 200 groups alarm events reviewing
- 4) 12 hours wave forms reviewing
- 5) Supporting for arrhythmia analysis and drug dose calculation

Standard:

Parameters: ECG, NIBP, SpO₂, RESP, TEMP, PR

Accessories: ECG cable, NIBP tube and cuff, SpO₂ Sensor, Temp probe, ECG Electrode, Ground cable, Power cable

Optional:

2-IBP, EtCO₂, Touch screen, Build in thermal printer