

8.4" Animal Monitor Model: Moni 8HV



Feature:

- 1) 8.4" color TFT LCD screen
- 2) Standard Configuration: Main unit with standard 6-parameter veterinary accessories (ECG, NIBP, SpO₂, RESP, TEMP, HR/PR), Rechargeable battery.
- 3) Optional Configuration: 2-IBP, EtCO₂ (Mainstream or Sidestream), Built-in Printer, Mobile Trolley, Wall Mount, 12V Power Supply.
- 4) Multi-language versions: English, Spanish, French, Russian, Arabic, Portuguese, Italian, Turkish, Chinese, etc.

10.4" Animal Monitor Model: Moni 8CV



Feature:

- 1) 10.4" color TFT LCD screen
- 2) Standard Configuration: Main unit with standard 6-parameter veterinary accessories (ECG, NIBP, SpO₂, RESP, TEMP, HR/PR), Rechargeable battery.
- 3) Optional Configuration: 2-IBP, EtCO₂ (Mainstream or Sidestream), CO, AG (Anesthesia Multigas Module), Built-in Printer, Mobile Trolley, Wall Mount, Touch Screen, 12V Power Supply for Ambulance Use, Carrying Bag, Wireless Network.
- 4) Multi-language versions: English, Spanish, French, Russian, Arabic, Portuguese, Italian, Turkish, Chinese, etc.

15" Animal Monitor Model: Moni 8LV



Feature:

- 1) 15" color TFT LCD screen
- 2) Remote controllable
- 3) Standard Configuration: Main unit with standard 6-parameter veterinary accessories (ECG, NIBP, SpO₂, RESP, TEMP, HR/PR), Rechargeable battery.
- 4) Optional Configuration: 2-IBP, EtCO₂ (Mainstream or Sidestream), CO, AG (Anesthesia Multigas Module), Built-in Printer, Mobile Trolley, Wall Mount, Touch Screen, 12V Power Supply for Ambulance Use, Carrying Bag, Wireless Network.
- 5) Multi-language versions: English, Spanish, French, Russian, Arabic, Portuguese, Italian, Turkish, Chinese, etc.

