


宠物笼

Pet Professional Cage

CCG-1



Cage Divider



Plastic Tray



Handle



Universal Wheel

- 1 可拆卸式门、隔板、网架、托盘使清理工作极其方便。  
Removable doors, Divider, Grid, Tray for easy Cleaning.
- 2 宠物笼材质为不锈钢使其持久耐用。(宠物笼门为“304”不锈钢, 侧板为“202”不锈钢材质)  
It's made by stainless steel for durability and long life-span. (Its door is made by "304" stainless steel, surrounds plate is made by "202" stainless steel).
- 3 宠物笼底部配备万向轮使其移动方便快捷。  
It's equipped with the universal wheel which it's easy to be moved.
- 4 小号笼尺寸 (CCG-1S): 长92\*深56\*高50cm      Small Cage Dimensions (CCG-1S): L92\*D56\*H50cm  
小号笼尺寸 (CCG-1XS): 长46\*深56\*高50cm      Small Cage Dimensions (CCG-1XS): L46\*D56\*H50cm  
中号笼尺寸 (CCG-1M): 长115\*深70\*高70cm      Medium Cage Dimensions (CCG-1M): L115\*D70\*H70cm  
大号笼尺寸 (CCG-1L): 长115\*深70\*高80cm      Large Cage Dimensions (CCG-1L): L115\*D70\*H80cm

[www.shivet.com](http://www.shivet.com)