

## Dual Color Lights Dental Curing Light

## Model: LC-C240C

## Features:

- 1) 5W big power light
- 2) 2 color lights: blue light for resin material solidification, white light for diagnosis function (or purple light for decayed tooth detection)
- 3) Four working models: strong power, low power, flashing , gradually strong
- 4) Solidify time: 5s~40s select
- 5) Wireless , low battery alert , battery can be recharged
- 6) Full charging to work for 10 seconds each time can be used continuously for more than 300 times
- 7) The base matched with light meter
- 8) It can solidify 23mm deepness in a few seconds
- 9) Color: black, blue, pink color is available



## Specifications:

Working voltage : 100-240V, 50/60Hz  
Wavelength: 420~480nm  
Light power: more than 1600mw/cm<sup>2</sup>  
Capability of battery: 2200mA/3.7V  
Solidification time and depth: 5s/3mm

## Standard:

Config 1: blue light + white light  
Config 2: blue light + purple light



Decayed Teeth



Healthy Teeth

