

Electric Operating Table

Model: ET-1200A

Features:

- 1) Electric lift operating table
 - 2) The table is made of 304 stainless steel, and equipped with stainless steel grid.
 - 3) Pedal lift control Easy to operate, stable and smooth.
- Size: 1200*600*(530-1000)mm



Electric Operating Table

Model: ET-1300A



Features:

Electric stainless steel X-type operating table with net is safe and reliable for users to operate.

It is made of high quality 304 stainless steel counter-tops.

With the bracing style elevator, it moves smoothly and steadily, and it's locked securely.

It is easy to operate with the touch-lift switch.

Size: 1300*600*(400-1100)mm



Electric Operating Table

Model: ET-1400A

Features:

- 1) The tabletop lifting is controlled by the electric foot pedal.
- 2) The working tabletop can be tilted 0°-15° left and right by the mechanical manual control.
- 3) The working tabletop can be tilted 0°-45° ahead and back by the mechanical manual control.
- 4) The tabletop and base is fully made of stainless steel.
- 5) The tabletop temperature can be adjusted from indoor temperature to 50° freely.
- 6) The machine has compact structure, reliable performance, reasonable manipulation convenient.



Specifications:

Tabletop size: 1400mm x 650mm

Tabletop height: 760-1060mm

Tabletop tilt angle: $\pm 15^\circ$

Joint space (stable and reliable): 0°