

Halogen Series

Model: LP02-5, LP01-5



LP02-5



LP01-5

Specifications:

Model	Illuminance	Colour temperature	Colour reduction index(Ra)	Illuminance depth	Total irradiance	Size of light field	Lamp head diameter	(24V) Rated power of bulb	Power supply voltage	Lowest height of installation	Bulb
LP02-5	≥50,000 Lux	3500±500K	88	600mm	490 W/m ²	160-220mm	500mm	25W	AC 220V/50-60 Hz (110V/60Hz custom-made)	2900mm	Cold light source

Model	Illuminance	Colour temperature	Colour reduction index(Ra)	Illuminance depth	Total irradiance	Size of light field	Power supply voltage	(24V) Rated power of bulb	Bulb
LP01-5	≥60,000 Lux	3700±500K	85	900mm	547W/m ²	160-220mm	AC 220V/50-60 Hz (110V/60Hz custom-made)	25W	Cold light source