

Halogen Series

Model: LP200, LP200E, LP200C, LP200W



Specifications:

Model	Illuminance	Colour temperature	Power supply voltage	(24V) Rated power of bulb	Lamp head diameter	Bulb	Emergency time	Lowest height of installation
LP200	≥50,000 Lux	4000±500K	AC 220V/50-60 Hz (110V/60Hz custom-made)	50W	180mm	Tungsten halogen bulb	/	/
LP200E	≥50,000 Lux	4000±500K	AC 220V/50-60 Hz (110V/60Hz custom-made)	50W	180mm	Tungsten halogen bulb	≥2 hours	/
LP200C	≥50,000 Lux	4000±500K	AC 220V/50-60 Hz (110V/60Hz custom-made)	50W	180mm	Tungsten halogen bulb	/	2300mm
LP200W	≥50,000/Lux	4000±500K	AC 220V/50-60 Hz (110V/60Hz custom-made)	50W	180mm	Tungsten halogen bulb	/	2300mm