

Hydraulic Operating Table

Model: OT-01AV(u)

Features:

- 1) Noncorrosive table-board is made of stainless steel plate with high quality, which makes easy washing and disinfection.
- 2) Right-and-left gradient is 0-15°, regulated by the manual machinery system.
- 3) Base plate is made of A3 armor plate with baking finish, and the casters add its mobility.
- 4) Matched with the operating bracket tray.
- 5) Designed with reasonable machinery, credible performance and simple operation.
- 6) Table-board height is 830mm.
- 7) Table-board size is 1000mm*500mm.



Hydraulic Operating Table

Model: OT-02CV

Features:

- 1) Using the hydraulic control on the operating table surface level movement.
- 2) Right-and-left gradient is 15° respectively by mechanical manual operation.
- 3) Fore-and-aft gradient is 45° by manual manipulation.
- 4) Overall structure is compact, reasonably reliable, easy manipulation.
- 5) V-shaped design and a front-end downcomer, easy for cleaning and disinfection.
- 6) Make the table-board up or down by stepping on the pedal.
- 7) Base plate is made of the A3 armor plate with baking finish and the high-class Italian simulation skin.
- 8) Surgical microscopes and operating bracket tray can be installed.
- 9) Table-board height is from 820mm to 930mm.
- 11) Table-board size is 1400×650mm.



Hydraulic Operating Table For Animal Heart

Model: OT-02HV

Features:

- 1) Table designed in U-top with drain, hydraulic System to adjust table height.
- 2) 0°-15° at right-and-left side inclined with manual operation.
- 3) Top made by Stainless Steel material helps to prevent from corrosion, rust and high temperature.
- 4) Equipped with a hole for ultrasound scanner at the upper-left of the top (The operators can open or close the hole according to their demands).
- 5) Base made by A3 armor plate, and with wheels to move conveniently.
- 6) Table Size(L*W): 1400mm*600mm
- 7) Table Height: 820-930mm

