

Infrared Therapeutic Lamp

Model: IR-I, IR-II

Applications:

Infrared therapy lamp is applied to infrared therapy treatment in medical units or pet hospitals.

Specifications:

Power supply: 110-220V/50-60Hz

Power: 400W

Adjustable angle of lamp holder(for IR-I): 60°

Adjustable angle of lamp holder(for IR-II): 360°

Adjustable height of lamp holder(for IR-II): 50-135cm

Time controllable: mechanical timing: 0-60 minutes

electronic timing: 0-95 minutes (optional)

IR wavelength: 2-10 μ m



IR-I



IR-II