

## LED Series

## Model: LEDLP01-4, LEDLP01-5



LEDLP01-5



LEDLP01-4

## Specifications:

Model	Illuminance	Colour temperature	Colour reduction index(Ra)	Illuminance depth	Total irradiance	Size of light field	Power supply voltage	(24V) Rated power of bulb	Bulb
LEDLP01-5	≥80,000 Lux	5000±500K	92	500mm	275W/m <sup>2</sup>	160-220mm	AC 220V/50-60 Hz (110V/60Hz custom-made)	3.2V/1W	LED Bulb ≥30,000h
LEDLP01-4	≥70,000 Lux	5000±500K	92	500mm	275W/m <sup>2</sup>	160-220mm	AC 220V/50-60 Hz (110V/60Hz custom-made)	3.2V/1W	LED Bulb ≥30,000h