

LED Series

Model: LEDLP01-4E, LEDLP01-5E



LEDLP01-5E



LEDLP01-4E

Specifications:

Model	Illuminance	Colour temperature	Colour reduction index(Ra)	Illuminance depth	Total irradiance	Size of light field	Power supply voltage	(24V) Rated power of bulb	Bulb	Emergency time
LEDLP 01-5E	≥80,000 Lux	5000±500K	92	500mm	275W/m ²	160-220mm	AC 220V/50-60 Hz (110V/60Hz custom-made)	3.2V/1W	LED Bulb ≥30,000h	≥6 hours
LEDLP 01-4E	≥70,000 Lux	5000±500K	92	500mm	275W/m ²	160-220mm	AC 220V/50-60 Hz (110V/60Hz custom-made)	3.2V/1W	LED Bulb ≥30,000h	≥6 hours