

LED Series

Model: LEDLP01A, LEDLP01W



LEDLP01A



LEDLP01W

Specifications:

Model	Illuminance	Colour temperature	Size of light field	Brightness adjustment	Lamp head diameter	Power supply voltage	Rated power of bulb	Bulb life	Pole
LEDLP01A	30,000Lux 0.5m /10,000Lux 1m	4800±200K	100mm	Absolute adjustment	55mm	AC 110/220V 50-60Hz	3. 3V/3W	≥20,000 hours	Adjust up and down
LEDLP01W	30,000Lux 0.5m /10,000Lux 1m	4800±200K	100mm	Absolute adjustment	55mm	AC 110/220V 50-60Hz	3. 3V/3W	≥20,000 hours	/