



LEDLP02-5

## Specifications:

Model	Illuminance	Colour temperature	Colour reduction index(Ra)	Illuminance depth	Total irradiance	Size of light field	Lamp head diameter	Rated power of bulb	Bulb	Power supply voltage	Lowest height of installation
LEDLP02-5 On Wall	≥100,000LUX	4800±300K	92	≥500mm	480W/m <sup>2</sup>	160-220mm	500mm	3.2V/1W	LED Bulb ≥30,000h	AC110-240V, 50/60Hz	2300mm
LEDLP02-5+5 On Ceiling	≥100,000LUX ≥100,000LUX	4800±300K	92 92	≥500mm ≥500mm	480W/m <sup>2</sup> 480W/m <sup>2</sup>	160-220mm	500mm	3.2V/1W	LED Bulb ≥30,000h	AC110-240V, 50/60Hz	2600mm
LEDLP02-5 On Ceiling	≥100,000LUX	4800±300K	92	≥500mm	480W/m <sup>2</sup>	160-220mm	500mm	3.2V/1W	LED Bulb ≥30,000h	AC110-240V, 50/60Hz	2300mm