

LED Series

**Model: LEDLP200, LEDLP200E,
LEDLP200C, LEDLP200W**

LEDLP200



LEDLP200E



LEDLP200C



LEDLP200W



Specifications:

Model	Illuminance	Colour temperature	Power supply voltage	Rated power of bulb	Lamp head diameter	Emergency time	Lowest height of installation
LEDLP200	≥50,000 Lux	6000±700K	AC 220V/50-60 Hz (110V/60Hz custom-made)	3.2V/1W	180mm	/	/
LEDLP200E	≥50,000 Lux	6000±700K	AC 220V/50-60 Hz (110V/60Hz custom-made)	3.2V/1W	180mm	≥6 hours	/
LEDLP200C	≥50,000 Lux	6000±700K	AC 220V/50-60 Hz (110V/60Hz custom-made)	3.2V/1W	180mm	/	2300mm
LEDLP200W	≥50,000/ Lux	6000±700K	AC 220V/50-60 Hz (110V/60Hz custom-made)	3.2V/1W	180mm	/	2300mm