

Veterinary Monitor (Bluetooth)

<http://www.shivet.com>

Feature

-  One button, easy to operate.
-  Artful design, easy to take with.
-  Build in rechargeable battery.
-  Multi parameter real time monitoring.
-  TFT display wave line display data.
-  Transmitting data to all kinds of digital devices via Bluetooth.
-  Suitable for: pet hospital, zoo, animal hospital, pet family.

Model: Moni 2A (Vet)

-  Blood Pressure
-  SPO2
-  Sphygmus
-  ECG
-  Temperature
-  Bluetooth



 For Android

Patient Monitor



Brand	Shinova
Type	Moni 2A
Display	2.4 inch color TFT LCD
ECG	Waveform: Lead I Band width 0.05Hz-100Hz(3dB) diagnose mode 1-40Hz monitor mode 0.1-25Hz operation mode
Temp	Type: Vascular/surface temperature Scope: 20°C-50°C Accuracy: $\pm 0.2^{\circ}\text{C}$
NIBP	Detection scope: Systolic blood pressure: 4kPa-34kPa(30-255mmHg) Diastolic pressure: 2kPa-29.3kPa(15-220mmHg) Mean pressure: 2.7kPa-31.3kPa (20-235mmHg) Detection accuracy: Static pressure: $\pm 4\text{mmHg}$
SpO2	Detection scope: 35%-100% Accuracy: $\pm 2\%$ 80%-100% $\pm 3\%$ 70%-79%
Heart rate	Scope: 15bpm to 300bpm Accuracy: $\pm 1\text{bpm}$

