


Portable Dental Unit  
Model: PD-852



67

**Specifications:**

Power: 600W  
 Maximum pressure: 6.5 bar  
 Clean water bottle: 1 L  
 Drain bottle for strong suction : 1 L  
 Drain bottle for weak suction : 0.5 L  
 Air bottle: 1 L  
 G.W./N.W: 34/29 KG  
 Dimension: 550\*420\*760 mm

**Standard:**

4 hole or 2 hole handpiece tube 2pcs  
 3-way syringe 1pc  
 Saliva ejector 1pc  
 High volume vacuum 1pc  
 LED curing light 1pc  
 Ultrasonic Scaler 1pc  
 Compressor head (550W) 1 set  
 On/Off Foot Switch 1pc  
 Water Supply 1 set

**Optional:**

Handpieces 2 pcs  
 Fiber optic handpiece 2 pcs  
 Mobile patient chair 1 pc  
 Mobile dentist stool 1 pc  
 Operating light 1 pc