

Portable Dental Unit Model: PD-892



Specifications:

Power: 550W
 Maximum pressure: 6.5 bar
 Clean water bottle: 1 L
 Drain bottle: 1 L
 Air aluminium tank: 6 L
 G.W./N.W: 33/28 KG
 Dimension: 550*420*760 mm

Standard:

4 hole or 2 hole handpiece tube 2 pcs
 3-way syringe 1 pc
 Saliva ejector 1 pc
 Compressor head (550W) 1 set
 On/Off Foot Switch 1 pc
 Water Supply 1 set

Optional:

Handpieces 2 pcs
 Fiber optic handpiece 2 pcs
 Mobile patient chair 1 pc
 Mobile dentist stool 1 pc
 Operating light 1 pc