

Portable Dental Unit
Model: PD-893

**Specifications:**

Power: 600W
Maximum pressure: 6.5 bar
Clean water bottle: 1 L
Drain bottle: 1 L
Air aluminous tank: 6 L
G.W./N.W: 34/29 KG
Dimension: 550*420*760 mm

Standard:

4 hole or 2 hole handpiece tube 2 pcs
3-way syringe 1 pc
Saliva ejector 1 pc
LED curing light 1 pc
Ultrasonic Scaler 1 pc
Compressor head (550W) 1 set
On/Off Foot Switch 1 pc
Water Supply 1 set

Optional:

Handpieces 2 pcs
Fiber optic handpiece 2 pcs
Mobile patient chair 1 pc
Mobile dentist stool 1 pc
Operating light 1 pc