

## Portable Veterinary Anesthesia Machine

### Model: AneCompact

#### Features:

AneCompact is a compact veterinary anesthesia machine that is very stable and reliable. It is widely used for pet hospital, veterinary clinic, animal zoo, farm department and so on.

#### Specification:

Respiration mode: Manual mode, open mode, closed mode, semi-open and semi-closed mode  
 Gas source: O<sub>2</sub>, quick inlet port  
 O<sub>2</sub> FLUSH: 20L/min~75L/min  
 Flow meter: O<sub>2</sub>, 0L/min ~10L/min, adjustable easily  
 Airway pressure gauge: -20cmH<sub>2</sub>O-100cmH<sub>2</sub>O  
 APL valve (Pop Off valve): adjustable easily, range: 0-60cmH<sub>2</sub>O  
 Inspiration valve and expiration valve: visible dome-type, ISO standard port, diameter 22mm  
 Soda lime chamber: 2 Liter, easy to change soda lime within one minute  
 Bains connector: open breath system, like Bains circuit and Jackson circuit  
 Anesthesia vaporizer: Enflurane, Isoflurane, Halothane or Sevoflurane vaporizer. Supply selectatec compatible mounting bar, ISO 23mm mounting bar is also compatible.  
 Auxiliary gas source output: supply the gas source to anesthesia ventilator  
 Three methods: tabletop, trolley or wall-mounted type (trolley and wall-mounted type is optional)



Wall-mounted type (optional)

#### Standard configuration:

Anesthesia machine, one vaporizer, one disposable breath circuit, one reservoir bag, one O<sub>2</sub> gas source pipeline.

#### Optional configuration:

Enflurane, Isoflurane, Halothane or Sevoflurane vaporizer; Trolley, wall-mounted arm, shelf for animal monitor.



Stand type (optional)

## Anesthesia Machine For Large Animals

### Model: AneKing

#### Features:

AneKing is a professional large animal anesthesia machine designed and developed by Shinova. It is specially used for large animal, like horse, cow, tiger, pig etc. It can be used in farm, zoo, wild animal zoo, forest park, animal research center and so on.

#### Specification:

Inspiration mode: open mode, closed mode, semi-open and semi-closed mode  
 Gas source: O<sub>2</sub>, 0.3~0.4Mpa  
 O<sub>2</sub> Flush valve: 25L/min~75L/min  
 Flow meter: high precision O<sub>2</sub> flow meter, 0~1Lmin~10Lmin  
 CO<sub>2</sub> circle absorber: with soda lime chamber, APL valve, airway pressure gauge, inspiration valve and expiration valve, the inspiration and expiration valve is demo-type and visible  
 APL valve: 0~60cmH<sub>2</sub>O, adjust easily  
 Airway pressure gauge: -20cmH<sub>2</sub>O~100cmH<sub>2</sub>O, -2Kpa~10Kpa  
 Soda lime chamber: 8 Liter, it is easy to change soda lime and used for over 1000kg animal  
 Bains connector: it is used for open breath circuit, connect with Jackson breath circuit or Bains breath circuit  
 Anesthesia vaporizer: Enflurane, Isoflurane, Halothane, Sevoflurane vaporizer with selectatec compatible mounting bar.  
 Cage mounting vaporizer is compatible  
 A handle on the top  
 Trolley is easily to transport everywhere (optional)



Trolley (optional)

#### Standard configuration:

Anesthesia machine, one vaporizer, one gas source pipeline, one set disposable breath circuit, one airbag.

#### Optional configuration:

O<sub>2</sub>&N<sub>2</sub>O two tubes flow meter, Trolley.