

Portable Veterinary Color Doppler Model: DopScan 6100V



Features:

- 1) Monitor supporting 180 °free adjustment, a wide-ranged size choice, 19 inch LCD monitor standard
- 2) Operation panel characterizing numerical control back light, fast regulation and one-key easy optimization
- 3) Four probe connectors, free probe switching, hot-swapping available
- 4) DVD-RW, automatically starting U disk with system restoration and maintenance function
- 5) Newly designed ergonomic console
- 6) Multiple choices of language display
- 7) DICOM3.0, remote network assistance available
- 8) Graphic workstations, import from and export to typical cases generating quick diagnosis

Clinical Application:

For clinical diagnosis of horses, cattle, sheep, pigs, cats, dogs and other animals.

Standard:

Main unit + 6.5MHz Rectal Linear Probe + 3.5Mhz
Convex probe;
15 inch LCD monitor, 3-probe connectors, Inner
Rechargeable Battery

Optional:

Probe: Vet Endo-rectal probe, Micro-convex Probe,
Convex Probe, Linear Probe, Phased Probe, Volume
Probe
Biopsy Bracket
3D/4D Module
Thermal Printer
Foot Switch
Trolley

