


Portable Veterinary ECG Machine

Model: ECG-T3V

Feature:

- 1) 3 or 1 channel ECG lead display
- 2) 6.2" smart touch screen with slide and multi-finger operate function
- 3) 122 kinds of automatic arrhythmia detection and printing
- 4) 40,000 ECG built-in records with binary image format
- 5) Real-time sampling and 600 seconds waveforms review
- 6) Analysis & print for any ECG leads within 10 seconds
- 7) Automatic standby and timing on-off function
- 8) Unlimited USB data transfer
- 9) Smart fuzzy ID search function
- 10) Ultra-quick mode battery recharge
- 11) Built-in thermal printer
- 12) Built-in Li-battery for up to 7 hours durability



Standard:

- Veterinary 5-Lead Clip ECG Cable
- Recording Paper (roll, 80mm*20m)
- Paper Roller
- Power Cord
- Fuse
- Ground Cable
- User Manual

Optional:

- Trolley
- Recording Paper
- External printer
- Built-in WiFi
- Scanner connectable
- Connectable with FTP server for data remote-control
- Connectable with PC-based ECG management software

Handheld Veterinary Pulse Oximeter

Model: OXI-90V

Feature:

- 1) Digital Technology
- 2) Alarm with freely set limit.
- 3) Accurate measurement during low perfusion
- 4) Anti-Movement, Stable Performance
- 5) Built-in rechargeable battery and seat-type charger; Power-Saved Function
- 6) Built-in memory card for data memory; connection with computer for downloading data
- 7) Bright, Easy-to-Read LED Numerals
- 8) Lightweight, Portable
- 9) 2 parameters: SpO₂ + Pulse Rate

Specifications:

SpO₂
 Measurement range: 0-100%
 Accuracy: ±2% during 70%-100%
 ±3% during 60%-69%

Pulse Rate

Measurement range: 15bpm-400bpm%
 Accuracy: 1 bpm or ±2%, Take the bigger one

Visual Displays and Indicators:

SpO₂ display: Red LED, Variable brightness
 Pulse Rate display: Green LED, variable brightness
 Pulse intensity: Red LED 8 bars
 Power Indicator Light:
 Normal: Green
 Low Power voltage: Red and Blink
 Low perfusion indicator: Red bar

Audio Indicator:

Probe off audio indicator (every thirty seconds)
 Pulse audio indicator
 Key-pressing indicator
 Low power voltage

Environmental Specifications:

Temperature Range: Operating: (0 ~ +50) °C
 Transport/Storage: (-40 ~ +70) °C
 Relative Humidity: Operating: 10~ 95%
 (noncondensing) Transport/Storage: 0~95%

Power Requirements:

Input: 98~260 VAC, 50/60Hz
 Consumption: < 80VA

Keys:

Pellicle keys

Standard Configuration:

Main unit, veterinary SpO₂ probe, Wrist strap, Rechargeable battery, Charger Pedestal, Power cord, Communication cable, CD-Rom, Operation manual

