

Tabletop Steam Sterilizer

Model: MS-TB20J, MS-TB24J

Features:

- 1) Indicator light indicates working state.
- 2) For 4~6 minutes rapidly sterilizing.
- 3) Sterilizing temperature and time can be preset.
- 4) Steam-water inner circulation system: no steam discharge, and the environment for sterilizing will be clean and dry.
- 5) Over temperature & over pressure auto-protection.
- 6) Safe protection of water lacking.
- 7) Exhaust the cool air automatically.
- 8) Door safety lock system.
- 9) With three stainless steel sterilizing plates.
- 10) Fully stainless steel structure.
- 11) Automatically shut off with beep reminding after sterilization.



Specification:

Model	MS-TB20J	MS-TB24J
Sterilizing chamber volume	20L (φ250×420mm)	24L (φ250×520mm)
Maximum working pressure	0.22Mpa	
Maximum working temperature	134 °C	
Adjustment of temperature	105-134°C	
Timer	0-60 min	
Chamber temperature equal	≤±1°C	
Source power	1.5KW, AC 220V/50Hz (or AC 110V/60Hz)	
Sterilizing plate	340×200×30mm (3 pieces)	400×200×30mm (3 pieces)
Dimension	480×480×384mm	580×480×384mm
Package dimension	700×580×500mm	800×580×500mm
G. W / N.W	43/38kg	50/44kg