


Ultrasonic Scaler
Model: US-A8



Specifications:

Color: white
 Vibration frequency of the tip: 25~31KHz
 Primary vibration excursion of the tip: $\leq 100\mu\text{m}$
 Output half-excursion: $< 2\text{N}$
 Ultrasonic output power: 3W~20W
 Water pressure: 0.01MPa~0.05MPa
 Input: AC100~240V, 47~63Hz
 Main unit input: DC30V 1.5A
 Foot switch battery: 1.5V \times 2
 Receiving sensitivity: -114Db
 Power : 30VA~48VA
 Weight of main unit: 1.10kg
 Weight of adapter: 0.25kg
 Dimensions: 330mm \times 155mm \times 155mm
 Gross Weight: 3.50kg
 Packing Size: 340 \times 240 \times 180mm

Features:

Luminous Detachable Handpiece
 Wireless control foot switch
 Scaling & Endo & Perio, Automatic Water Supply