

## Veterinary Bone Drill

### Model: DS-0110

#### Features:

- 1) Imported material, manufactured with advanced technology
- 2) Lower weight and noise, less operation time
- 3) Pakaged assembly and modified structure are easier to use
- 4) NI-MH battery without memory environment
- 5) Autoclavable



#### Specifications:

Hand piece: Working voltage: 14.4V; Output power:  $\geq 90W$ ;  
Unit rise in temperature:  $\leq 20^{\circ}C$   
Charger: Input voltage: AC 110V-240V, 50Hz-60Hz; Output  
voltage: 14.4V, 900mA; Input power: 30W  
Battery: Voltage power: 14.4V; Charging time: about 3  
hours  
Speed: 0-1300 rpm  
Noise:  $\leq 65db$   
Sterilized:  $\leq 135^{\circ}C$

#### Standard Configuration:

1\*Hand piece, 2\*Battery, 1\*Charger, 1\*Drill chuck, 1\*Sterilizing channel, 1\*Key, 1\*Aluminium box.

## Veterinary Canulated Drill

### Model: DS-0210

#### Features:

- 1) Axial run-out of concentricity is within 0.05mm.
- 2) Rotating speed of 600 rpm is more suited to the medical standard for the operation power.
- 3) About 4 N/M for the foreign brand, and 3.3 N/M for Shinova product.
- 4) Supporting the interface of the brand product.
- 5) Full-sealed manufacturing technology for electric motor more complies with the requirements of high temperature and high pressure disinfection.
- 6) Semi-luster surface treatment technology gives the operator more comfortable operation.
- 7) Multipurpose interface is easier to use in clinic.
- 8) 4.0mm hollow for all brand guide needles to pass through.

#### Specifications:

Hand piece: Working voltage: 14.4V; Output power:  $\geq 90W$ ;  
Unit rise in temp:  $\leq 20^{\circ}C$   
Charger: Input voltage: AC 110V-240V, 50Hz-60Hz; Output  
power: 14.4V, 900mA; Input power: 30W  
Battery: Voltage: 14.4V; Charging time: 2 hours or more  
Speed: 0-650 rpm  
Noise:  $\leq 65db$   
Sterilized:  $\leq 135^{\circ}C$



#### Standard Configuration:

1\*Hand piece, 1\*Charger, 1\*Battery, 1\*Drill chuck, 1\*Sterilized channel, 1\*Key, 1\*Aluminum box.

## Veterinary Oscillating Saw

### Model: DS-0310

#### Features:

- 1) Long working time
- 2) Autoclavable, can be sterilized up to  $135^{\circ}C$
- 3) Adopt total newly electric machinery
- 4) NI-MH battery without memory environment
- 5) For joint surgery

#### Specifications:

Hand piece: Working voltage: 14.4V; Output power:  $\geq 90W$ ; Unit rise in temperature:  $\leq 20^{\circ}C$   
Charger: Input voltage: AC 110V-240V, 50Hz-60Hz; Output  
voltage: 14.4V, 900mA; Input power: 30W  
Battery: Voltage power: 14.4V; Charging time: 2 hours or  
more  
Speed: 0-20000 osc/m  
Noise:  $\leq 70db$   
Sterilized:  $\leq 135^{\circ}C$

#### Standard Configuration:

1\*Hand piece, 1\*Battery, 1\*Charger, 2\*Blade saw, 1\*Sterilizing channel, 1\*Aluminium box.

