

## Veterinary Fully-automatic Chemistry Analyzer

## Model: Chemo 120V

## Features:

- 1) Fully automatic biochemical analyzer with specialized veterinary software
- 2) 8 pre-determined animals settings and unlimited number of user-defined animal settings
- 3) Throughput up to 100 tests per hour
- 4) Labor saving, simple programming and real walk-away operation
- 5) Reagent open system, closed system on request
- 6) Micro-volume for sample and reagent
- 7) Sample position and reagent position switchable
- 8) Multi-language software supported

## Specifications:

Animals types: Dog, Cat, Cow, Fish, Horse, Mouse, Pig, Rabbit, etc.

Contents of testing: Clinic chemistry, Immunoturbidimetric

Analytical modes: End-point, Two Point, Kinetic

Programming: Open system with user defined programs

Automatic, Random Access

STAT sample priority

9 postions for sample

Sample Volume: 3~45uL, 0.1uL adjustable

26 postions for reagents

Refrigerated temperature: 4-15°C

Reagent volume: R1: 180-450uL, 0.1uL adjustable

R2: 3-250uL, 0.1uL adjustable

Reaction Postion: 90 cuvettes, 10 cuvettes/strip

Cuvette: optical length 5mm

Reaction volume: 180-500uL

Reaction temperature: 37°C

Lamp: Halogen lamp

Absorbance range: 0-3.500Abs

Wavelengths: 8 wavelengths from 340-700nm

Resolution: 0.0001Abs

Power supply: AC 110V/220V, 50Hz/60Hz

Water consumption: <3L/hour

Dimension L\*W\*H(mm): 630\*480\*510

Net weight: 70kg

