

Veterinary Micro Bone Drill

Model: DS-1410

Features:

- 1) The global leading technology, advanced art of manufacturing, the smallest hand-piece and precision.
- 2) All of these advantages of the product can meet various requirements of orthopedic surgery and microsurgery.
- 3) Core drilling speed of 1100 rpm.
- 4) Swing saw speed of 15000 times per min.

**Specifications:**

Hand Piece: Working Voltage: 7.2V; Output Power: $\geq 50W$; Unit rise in temperature: $\leq 20^{\circ}C$
Charger: Input Voltage: AC 110V-240V, 50Hz-60Hz; Output Power: 7.2V; Input Power: 30W
Battery: Voltage: 7.2V; Charging time: About 3 hours
Drill: Frequency: 1100 r/m; Non-load noise: $\leq 45dB$; Output power: 50W; Sterilized: $\leq 135^{\circ}C$

Standard Configuration:

1*Hand piece, 2*Battery, 1*Charger, 1*Drill chuck, 2*Sterilizing channel, 1*Aluminium case.