

## Veterinary Micro/Mini Multifunctional Tool

### Model: DS-0610A

#### Features:

- 1) Quick chuck: ensures outstanding performance during operation.
- 2) Oscillating saw: 360° Rotating head can be applied to any surgical Sites. Wrenchless collect allows for easy, fast and safe loading and unloading blades.
- 3) Mainframe is made of Aluminium, light weight, durable, long working time
- 4) Perfect ergonomic design.
- 5) Can meet different kinds of surgery operation.

#### Specifications:

Hand piece: Working Voltage: 7.2V; Output Power: ≥ 50W; Unit rise in temperature: ≤20°C;  
 Charger: Input Voltage: AC 110V-240V,50Hz-60Hz;  
 Output Power: 7.2V; Input Power: 30W;  
 Battery: Voltage: 7.2V; Charging time: About 2 hours or more;  
 Quick chuck (bone drill speed) attachment: Frequency: 0-1300 r/m; Noise: ≤60db; Torque: ≥1.6NM; Sterilize: ≤ 135°C;  
 Oscillating saw attachment: Frequency: 0-20000 osc/m; Noise: ≤70db; Sterilize: ≤135°C.

#### Standard Configuration:

1\*Hand piece; 1\*Charger; 2\*Battery 7.2V; 1\*Drill chuck;  
 2\*Sterilized channel; 2\*Saw blade; 2\*Key; 1\*Aluminum box; 1\*Quick chuck (bone drill speed) attachment;  
 1\*Oscillating saw attachment.



## Medical Electric Plaster Saw

### Model: PS-5001

#### Features:

- 1) The basic surgical equipment which can effectively cut all kinds of bandages.
- 2) High-speed saw blade vibration system.

#### Specifications:

- 1) Non-load speed: 12500r/min
- 2) Input power: 115VA
- 3) Noise: ≤84dB
- 4) Input voltage: 110-240V, 50Hz/60Hz

#### Standard Configuration:

1\*Saw unit, 3\*Saw blades, 1\*Allen wrench to load & unload saw blades, 1\*Instruction manual, 1\*Aluminum box.

