

Veterinary Operating Table for Large Animals

Model: OT-03BV

Features:

1. The operating table is powered by electric controlling and hydraulic lifting system.
2. The whole structure is compact and reliable, and it's very easy to operate.
3. The tabletop is fully made of stainless steel, which makes solid, anticorrosive, rust-proof and high-temperature resistance.
4. Stainless steel infusion pole is included.
5. The tabletop can be tilted: 0-30°C.
6. It can support 1500 kgs.
7. Equipped with universal brake wheel, it can move everywhere easily.
8. Table length and width: 3700mm(L) * 1200mm(W) (including pillow length 700mm)
9. Table height range from the ground: 600mm-870mm
10. Length and width of table pedestal: 700mm(L) * 400mm(W)

Optional:

- a. Heating System
- b. Customized size is available



