

Veterinary Ultrasonic Cleaner

Model: SB-50, SB-80, SB-100D, SB-120D, SB-3200D, SB-5200D, SB-4200D, SB25-12D

Features:

- 1) Digital display and setting ultrasonic cleaning time
- 2) Display working countdown time
- 3) Timing range: 1-600 minutes adjustable
- 4) Equipped with a special stainless steel basket, reduce sound cover
- 5) The shell and cover are made by stainless steel
- 6) The equipment controlled by Single-chip controller
- 7) The last working parameters can be stored in memory after power-off



SB25-12D



SB-50



SB-120D



SB-5200D

Specification:

Model	Inner Tank L*W*H (mm)	Volume (L)	Frequency (KHz)	Power (W)	Drainage	Adjustable Time (min.)	N. W. (Kg)	Packing Size (mm)	Rack
SB-50	150*85*65	0.8	40	50	/	1-120	2	220*150*210	√
SB-80	155*140*100	2	40	80	/	1-120	2.2	240*230*300	√
SB-100D	300*150*100	4.5	40	100	/	1-600	4.5	400*240*435	√
SB-120D(3L)	240*140*100	3	40	120	/	1-600	4.5	325*230*300	√
SB-120D(5L)	240*140*150	5	40	120	/	1-600	4	325*230*350	√
SB-3200D	300*155*150	6	40	180	/	1-600	5	400*240*380	√
SB-5200D	300*240*150	10	40	300	√	1-600	7.2	450*325*375	√
SB-4200D	300*240*200	14.4	40	400	√	1-600	9	445*320*460	√
SB25-12D	500*300*150	22.5	40	600	√	1-600	13.5	640*395*375	√